

Stimulating Electrodes Featuring Silver Technology



Our Advantrode® Elite electrodes feature our uniquely-applied layer of pure silver, optimal to provide the most comfortable and effective stimulation. Advantrode Elite electrodes are offered in a variety of flexible, durable topcoat materials: White Foam, Tan Tricot, Tan Spunlace and White Spunlace. PMI's state-of-the-art manufacturing processes, superior medical-grade materials and commitment to quality is reflected in each electrode. The expanded selections in different sizes, shapes and hydrogels provide you and your customers the freedom to choose the most effective electrodes for maximum treatment results.

Feature	Benefit
Highly conductive silver-coated carbon film	Provides uniform electrical current flow over the entire surface of the electrode.
35-mil dual layer hydrogel	Maximum reusability throughout the life of the electrode.
Tamper-evident bag *(Available with most item numbers)	Assurance that the product has not been altered.
Tin-coated copper wire with molded grip	Enhanced current dispersion with ultra-low impedance
Soft, molded grip wire socket	Ease of connection and disconnection, while providing superior comfort.
Large selection of shapes and sizes	Adaptable to multiple electrode application sites.
Foam and cloth topcoats	Ultimate flexibility and conformability to body contours.
Skin-friendly blue hydrogel	Designed for sensitive/delicate skin



Scan the QR Code to find out more information on Advantrode® Elite products



AdvanTrode® Elite Tan Spunlace

Item #	Description	Elec/Pack	Pk/Case
WTS2	2" (5cm) Round	4	10
WTS15	1.5" x 1.5" (4cm x 4cm) Square	4	10
WTS22	2" x 2" (5cm x 5cm) Square	4	10
WTS13	1.75" x 3.75" (4.5cm x 9.5cm) Rectangle	4	10



AdvanTrode® Elite White Foam

SF375	1.375" (3.5cm) Round w/Snap	4	10
WF2	2" (5cm) Round	4	10
WF3	3" (7.5cm) Round	4	10
WF22	2" x 2" (5cm x 5cm) Square	4	10
WF125o	1.5" x 2.5" (4cm x 6.5cm) Oval	4	10
WF13	1.75" x 3.75" (4.5cm x 9.5cm) Rectangle	4	10
WF35	3" x 5" (7.5cm x 13cm) Rectangle	2	10



AdvanTrode® Elite White Spunlace

WWS2	2" (5cm) Round	4	10
WWS22	2" x 2" (5cm x 5cm) Square	4	10
WWS13	1.75" x 3.75" (4.5cm x 9.5cm) Rectangle	4	10



AdvanTrode® Elite Tan Tricot

WT375	1.375" (3.5cm) Round	4	10
WT2	2" (5cm) Round	4	10
WT3	3" (7.5cm) Round	4	10
WT22	2" x 2" (5cm x 5cm) Square	4	10
WT24o	2" x 4" (5cm x 10cm) Oval	4	10
WT13	1.75" x 3.75" (4.5cm x 9.5cm) Rectangle	4	10
WT35o	3" x 5" (7.5cm x 13cm) Rectangle	2	10
WT47	4.25" x 7" (11cm x 18cm)	1	10



AdvanTrode® Elite White Foam Back Therapy

WF46	4" x 6" (10cm x 15cm) Butterfly	1	10
WF1315	13" x 1.5" (33cm x 4cm) Backlap	1	10



AdvanTrode® Elite Skin-Friendly Blue Gel

BWT11	1.75" x 1.75" (4.5cm x 4.5cm) Square	4	10
BWT13	1.75" x 3.75" (4.5cm x 9.5cm) Rectangle	4	10
BST11	1.75" x 1.75" (4.5cm x 4.5cm) Square w/Snap	4	10
BST13	1.75" x 3.75" (4.5cm x 9.5cm) Rectangle w/Snap	4	10